

## 兔抗 ORC4 多克隆抗体

中文名称：兔抗 ORC4 多克隆抗体

英文名称：Anti-ORC4 rabbit polyclonal antibody

别名：origin recognition complex subunit 4; ORC4L; ORC4P

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ORC4

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The origin recognition complex (ORC) is a highly conserved six subunit protein complex essential for the initiation of the DNA replication in eukaryotic cells. Studies in yeast demonstrated that ORC binds specifically to origins of replication and serves as a platform for the assembly of additional initiation factors such as Cdc6 and Mcm proteins. This gene encodes a subunit of the ORC complex. Several alternatively spliced transcript variants, some of which encode the same protein, have been reported for this gene.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	ORC4
<b>Immunogen:</b>	Fusion protein of human ORC4
<b>Full name:</b>	origin recognition complex subunit 4
<b>Synonyms:</b>	ORC4L; ORC4P
<b>SwissProt:</b>	O43929
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human cervical cancer
<b>IHC Recommend dilution:</b>	40-200



