

兔抗 VAMP5 多克隆抗体

中文名称: 兔抗 VAMP5 多克隆抗体

英文名称: Anti-VAMP5 rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

抗 原: VAMP5

宿 主: Rabbit

反应种属: Human

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Synaptobrevins/VAMPs, syntaxins, and the 25-kD synaptosomal-associated protein are the main components of a protein complex involved in the docking and/or fusion of vesicles and cell membranes. The VAMP5 gene is a member of the vesicle-associated membrane protein (VAMP)/synaptobrevin family and the SNARE superfamily. This VAMP family member may participate in vesicle trafficking events that are associated with myogenesis.
Applications:	ELISA, IHC
Name of antibody:	VAMP5
Immunogen:	Full length fusion protein

Full name:	vesicle-associated membrane protein 5
SwissProt:	O95183
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human prostate cancer and Human thyroid cancer
IHC Recommend dilution:	25-100



