

兔抗 RTKN2 多克隆抗体

中文名称: 兔抗 RTKN2 多克隆抗体

英文名称: Anti-RTKN2 rabbit polyclonal antibody

别名: rhotekin 2; PLEKHK1; bA531F24.1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: RTKN2

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	May play an important role in lymphopoiesis.
Applications:	ELISA, IHC
Name of antibody:	RTKN2
Immunogen:	Fusion protein of human RTKN2
Full name:	rhotekin 2
Synonyms:	PLEKHK1; bA531F24.1
SwissProt:	Q8IZC4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human liver cancer
IHC Recommend dilution:	50-300

