

兔抗 ORM2 多克隆抗体

中文名称：兔抗 ORM2 多克隆抗体

英文名称： Anti-ORM2 rabbit polyclonal antibody

别名： AGP2; AGP-B; AGP-B'

相关类别： 一抗

抗原： ORM2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.
Applications:	ELISA, WB, IHC
Name of antibody:	ORM2
Immunogen:	Fusion protein of human ORM2
Full name:	orosomuroid 2
Synonyms :	AGP2; AGP-B; AGP-B'

SwissProt:	P19652
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human lung cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	24 kDa
WB Positive control:	Mouse liver tissue
WB Recommended dilution:	200-1000



