

## 兔抗 IRAK2 多克隆抗体

中文名称:	兔抗 IRAK2 多克隆抗体
英文名称:	Anti-IRAK2 rabbit polyclonal antibody
别名:	interleukin 1 receptor associated kinase 2; IRAK-2
抗原:	IRAK2
储存:	冷冻 (-20℃)
宿主:	Rabbit
相关类别:	一抗
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	IRAK2
<b>Immunogen:</b>	Fusion protein of human IRAK2
<b>Full name:</b>	interleukin 1 receptor associated kinase 2
<b>Synonyms:</b>	IRAK-2
<b>SwissProt:</b>	O43187

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200

