

兔抗 OSR2 多克隆抗体

中文名称: 兔抗 OSR2 多克隆抗体

英文名称: Anti-OSR2 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃) 避光

宿 主: Rabbit

抗 原: OSR2

反应种属: Human Mouse Rat

标 记 物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).
Applications:	WB
Name of antibody:	OSR2
Immunogen:	Synthesized peptide derived from Internal of human OSR2.
Full name:	odd-skipped related transcription factor 2
Synonyms :	0
SwissProt:	Q8N2R0
WB Predicted band size:	36 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

