

兔抗 OTUD6A 多克隆抗体

中文名称: 兔抗 OTUD6A 多克隆抗体

英文名称: Anti-OTUD6A rabbit polyclonal antibody

别名: OTU deubiquitinase 6A; DUBA2; HSHIN6

相关类别: 一抗

储存: 冷冻 (-20℃)

抗原: OTUD6A

宿主: Rabbit

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Deubiquitinating enzymes (DUBs; see MIM 603478) are proteases that specifically cleave ubiquitin (MIM 191339) linkages, negating the action of ubiquitin ligases. DUBA2 belongs to a DUB subfamily characterized by an ovarian tumor (OTU) domain.
Applications:	ELISA, WB
Name of antibody:	OTUD6A
Immunogen:	Synthetic peptide of human OTUD6A
Full name:	OTU deubiquitinase 6A
Synonyms:	DUBA2; HSHIN6
SwissProt:	Q7L8S5

ELISA Recommended dilution:	5000-10000
WB Predicted band size:	33 kDa
WB Positive control:	Human breast cancer tissue, TM4, 231 and Jurkat cell lysates
WB Recommended dilution:	500-2000

