

兔抗 OTX1 多克隆抗体

中文名称：兔抗 OTX1 多克隆抗体

英文名称： Anti-OTX1 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： OTX1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the bicoid sub-family of homeodomain-containing transcription factors. The encoded protein acts as a transcription factor and may play a role in brain and sensory organ development. A similar protein in mouse is required for proper brain and sensory organ development and can cause epilepsy. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	OTX1
Immunogen:	Fusion protein of human OTX1
Full name:	orthodenticle homeobox 1
SwissProt:	P32242
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human liver cancer;Human esophagus cancer
IHC Recommend dilution:	20-100



