

## 兔抗 OVAL 多克隆抗体

中文名称：兔抗 OVAL 多克隆抗体

英文名称： Anti-OVAL rabbit polyclonal antibody

别名： SERPINB14

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： OVAL

宿主： Rabbit

反应种属： Chicken

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Ovalbumin (abbreviated OVA) is the main protein found in egg white, making up 60-65% of the total protein. Ovalbumin displays sequence and three-dimensional homology to the serpin superfamily, but unlike most serpins it is not a serine protease inhibitor. It is presumed to be a storage protein.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	OVAL
<b>Immunogen:</b>	OVAL protein
<b>Full name:</b>	ovalbumin (SERPINB14)
<b>Synonyms :</b>	SERPINB14

<b>SwissProt:</b>	P01012
<b>ELISA Recommended dilution:</b>	2000-5000
<b>WB Predicted band size:</b>	43 kDa
<b>WB Positive control:</b>	OVAL protein
<b>WB Recommended dilution:</b>	500-2000

