

兔抗 OVOL1 多克隆抗体

中文名称：兔抗 OVOL1 多克隆抗体

英文名称： Anti-OVOL1 rabbit polyclonal antibody

别名： HOVO1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： OVOL1

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a putative zinc finger containing transcription factor that is highly similar to homologous protein in Drosophila and mouse. Based on known functions in these species, this protein is likely involved in hair formation and spermatogenesis in human as well.
Applications:	WB
Name of antibody:	OVOL1
Immunogen:	Synthesized peptide derived from N-terminal of human OVOL1.
Full name:	ovo-like zinc finger 1
Synonyms :	HOVO1

SwissProt:	O14753
WB Predicted band size:	30 kDa
WB Positive control:	K562 cells and A549 cells lysates
WB Recommended dilution:	500-3000

