

兔抗 OXR1 多克隆抗体

中文名称: 兔抗 OXR1 多克隆抗体

英文名称: Anti-OXR1 rabbit polyclonal antibody

别名: TLDC3; Nbla00307

抗原: OXR1

储存: 冷冻 (-20℃) 避光

宿主: Rabbit

反应种属: Human Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	May be involved in protection from oxidative damage.
Applications:	IHC
Name of antibody:	OXR1
Immunogen:	Synthesized peptide derived from internal of human OXR1.
Full name:	oxidation resistance 1
Synonyms :	TLDC3; Nbla00307
SwissProt:	Q8N573
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100

