

## 兔抗 OXSR1 多克隆抗体

中文名称: 兔抗 OXSR1 多克隆抗体

英文名称: Anti-OXSR1 rabbit polyclonal antibody

别 名: OSR1

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: OXSR1

反应种属: Human, Mouse

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

Background:	The product of this gene belongs to the Ser/Thr protein kinase family of proteins. It regulates downstream kinases in response to environmental stress, and may play a role in regulating the actin cytoskeleton.
Applications:	ELISA, WB
Name of antibody:	OXSR1
Immunogen:	Fusion protein of human OXSR1
Full name:	oxidative stress responsive 1
Synonyms:	OSR1
SwissProt:	O95747
<b>ELISA Recommended dilution:</b>	1000-2000



WB Predicted band size:	58 kDa
WB Positive control:	Human normal liver tissue
WB Recommended dilution:	200-1000

