

## 兔抗 P2RY1 多克隆抗体

中文名称：兔抗 P2RY1 多克隆抗体

英文名称： Anti-P2RY1 rabbit polyclonal antibody

别名： P2Y1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： P2RY1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The product of this gene belongs to the family of G-protein coupled receptors. This family has several receptor subtypes with different pharmacological selectivity, which overlaps in some cases, for various adenosine and uridine nucleotides. This receptor functions as a receptor for extracellular ATP and ADP. In platelets binding to ADP leads to mobilization of intracellular calcium ions via activation of phospholipase C, a change in platelet shape, and probably to platelet aggregation.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	P2RY1

<b>Immunogen:</b>	Synthetic peptide of human P2RY1
<b>Full name:</b>	purinergic receptor P2Y, G-protein coupled, 1
<b>Synonyms :</b>	P2Y1
<b>SwissProt:</b>	P47900
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human breast cancer and Human colon cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	A549 cells and mouse liver tissue
<b>WB Recommended dilution:</b>	200-1000



