

## 兔抗 PACSIN2 多克隆抗体

中文名称：兔抗 PACSIN2 多克隆抗体

英文名称： Anti-PACSIN2 rabbit polyclonal antibody

别名： SDPII

相关类别： 一抗

储存： 冷冻（-20℃）

抗原： PACSIN2

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the protein kinase C and casein kinase substrate in neurons family. The encoded protein is involved in linking the actin cytoskeleton with vesicle formation by regulating tubulin polymerization. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PACSIN2
<b>Immunogen:</b>	Fusion protein of human PACSIN2
<b>Full name:</b>	protein kinase C and casein kinase substrate in neurons 2
<b>Synonyms :</b>	SDPII
<b>SwissProt:</b>	Q9UNF0

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human lung cancer and human colon cancer
<b>IHC Recommend dilution:</b>	25-100

