

## 兔抗 PADI4 多克隆抗体

中文名称：兔抗 PADI4 多克隆抗体

英文名称： Anti-PADI4 rabbit polyclonal antibody

别名： PAD; PAD4; PDI4; PDI5; PADI5

抗原： PADI4

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of a gene family which encodes enzymes responsible for the conversion of arginine residues to citrulline residues. This gene may play a role in granulocyte and macrophage development leading to inflammation and immune response.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PADI4
<b>Immunogen:</b>	Fusion protein of human PADI4
<b>Full name:</b>	peptidyl arginine deiminase 4
<b>Synonyms:</b>	PAD; PAD4; PDI4; PDI5; PADI5
<b>SwissProt:</b>	Q9UM07

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human ovarian cancer
<b>IHC Recommend dilution:</b>	40-200
<b>WB Predicted band size:</b>	74 kDa
<b>WB Positive control:</b>	SKOV3 and 293T cell lysates
<b>WB Recommended dilution:</b>	500-2000



