

兔抗 PAFAH1B1 多克隆抗体

中文名称：兔抗 PAFAH1B1 多克隆抗体

英文名称： Anti-PAFAH1B1 rabbit polyclonal antibody

别名： MDS, LIS1, LIS2, MDCR, PAFAH

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PAFAH1B1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This locus was identified as encoding a gene that when mutated or lost caused the lissencephaly associated with Miller-Dieker lissencephaly syndrome. This gene encodes the non-catalytic alpha subunit of the intracellular Ib isoform of platelet-activating factor acetylhydrolase, a heterotrimeric enzyme that specifically catalyzes the removal of the acetyl group at the SN-2 position of platelet-activating factor (identified as 1-O-alkyl-2-acetyl-sn-glycerol-3-phosphorylcholine). Two other isoforms of intracellular platelet-activating factor acetylhydrolase exist: one composed of multiple subunits, the other, a single subunit. In addition,

	a single-subunit isoform of this enzyme is found in serum.
Applications:	ELISA, IHC
Name of antibody:	PAFAH1B1
Immunogen:	Synthetic peptide of human PAFAH1B1
Full name:	platelet-activating factor acetylhydrolase 1b, regulatory subunit 1 (45kDa)
Synonyms :	MDS, LIS1, LIS2, MDCR, PAFAH
SwissProt:	P43034
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human ovarian cancer and Human brain
IHC Recommend dilution:	25-100

