

## 兔抗 PAFAH1B2 多克隆抗体

中文名称：兔抗 PAFAH1B2 多克隆抗体

英文名称：Anti-PAFAH1B2 rabbit polyclonal antibody

别名：platelet activating factor acetylhydrolase 1b catalytic subunit 2; HEL-S-303

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PAFAH1B2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Platelet-activating factor acetylhydrolase (PAFAH) inactivates platelet-activating factor (PAF) into acetate and LYSO-PAF. This gene encodes the beta subunit of PAFAH, the other subunits are alpha and gamma. Multiple alternatively spliced transcript variants have been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PAFAH1B2
<b>Immunogen:</b>	Fusion protein of human PAFAH1B2
<b>Full name:</b>	platelet activating factor acetylhydrolase 1b catalytic subunit 2

<b>Synonyms:</b>	HEL-S-303
<b>SwissProt:</b>	P68402
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-300

