

## 兔抗 RTN1 多克隆抗体

中文名称：兔抗 RTN1 多克隆抗体

英文名称： Anti-RTN1 rabbit polyclonal antibody

别名： reticulon 1; NSP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RTN1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene belongs to the family of reticulon encoding genes. Reticulons are associated with the endoplasmic reticulum, and are involved in neuroendocrine secretion or in membrane trafficking in neuroendocrine cells. This gene is considered to be a specific marker for neurological diseases and cancer, and is a potential molecular target for therapy. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	RTN1
<b>Immunogen:</b>	Fusion protein of human RTN1

<b>Full name:</b>	reticulon 1
<b>Synonyms:</b>	NSP
<b>SwissProt:</b>	Q16799
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	84 kDa
<b>WB Positive control:</b>	Mouse brain tissue lysate
<b>WB Recommended dilution:</b>	500-2000



