

兔抗 RTN4R 多克隆抗体

中文名称：兔抗 RTN4R 多克隆抗体

英文名称：Anti-RTN4R rabbit polyclonal antibody

抗 原：RTN4R

别 名：reticulon 4 receptor NGR; NOGOR

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the receptor for reticulon 4, oligodendrocyte myelin glycoprotein and myelin-associated glycoprotein. This receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system.
Applications:	ELISA, IHC
Name of antibody:	RTN4R
Immunogen:	Fusion protein of human RTN4R
Full name:	reticulon 4 receptor
Synonyms:	NGR; NOGOR

SwissProt:	Q9BZR6
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human gastric cancer;

