

兔抗 PAK1(Ab-199) 多克隆抗体

中文名称: 兔抗 PAK1(Ab-199) 多克隆抗体

英文名称: Anti-PAK1(Ab-199) rabbit polyclonal antibody

别名: PAKalpha

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

抗原: PAK1(Ab-199)

宿主: Rabbit

反应种属: Human Mouse Rat

标记物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:	This gene encodes a family member of serine/threonine p21-activating kinases, known as PAK proteins. These proteins are critical effectors that link RhoGTPases to cytoskeleton reorganization and nuclear signaling, and they serve as targets for the small GTP binding proteins Cdc42 and Rac. This specific family member regulates cell motility and morphology. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	WB, IHC
Name of antibody:	PAK1(Ab-199)
Immunogen:	Synthesized non-phosphopeptide derived from human PAK

	1/2 around the phosphorylation site of serine 199 (T-K-S(p)-V-Y).
Full name:	p21 protein (Cdc42/Rac)-activated kinase 1
Synonyms :	PAKalpha
SwissProt:	Q13153
IHC positive control:	Human placenta tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	61 kDa
WB Positive control:	K562 cells lysate
WB Recommended dilution:	500-3000



