

兔抗 VASH2 多克隆抗体

中文名称: 兔抗 VASH2 多克隆抗体

英文名称: Anti-VASH2 rabbit polyclonal antibody

储存: 冷冻 (-20℃)

抗原: VASH2

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	Vasohibin-2, also known as VASH2, is a 355 amino acid protein belonging to the Vasohibin family. Expressed in various embryonic organs, Vasohibin-2 can be detected during embryonic weeks six through twelve and is found in vessels of 20-week embryonic organs as well as in endothelial cells of neonatal large vessels. Induced by VEGF, Vasohibin-2 is an angiogenesis inhibitor and interferes with the proliferation, migration and network formation by endothelial cells.
Applications:	ELISA, WB
Name of antibody:	VASH2
Immunogen:	Full length fusion protein

Full name:	vasohibin 2
SwissProt:	Q86V25
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	40 kDa
WB Positive control:	Hela and 293T cell
WB Recommended dilution:	500-2000

