

兔抗 RTP4 多克隆抗体

中文名称：兔抗 RTP4 多克隆抗体

英文名称： Anti-RTP4 rabbit polyclonal antibody

别名： receptor transporter protein 4; IFRG28; Z3CXXC4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： RTP4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Probable chaperone protein which facilitates trafficking and functional cell surface expression of some G-protein coupled receptors (GPCRs). Promotes functional expression of the bitter taste receptor TAS2R16 (PubMed:16720576). Also promotes functional expression of the opioid receptor heterodimer OPRD1-OPRM1 (By similarity).
Applications:	ELISA, IHC
Name of antibody:	RTP4
Immunogen:	Fusion protein of human RTP4
Full name:	receptor transporter protein 4

Synonyms:	IFRG28; Z3CXXC4
SwissProt:	Q96DX8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-200

