

## PLIN3 抗原（重组蛋白）

中文名称： PLIN3 抗原（重组蛋白）

英文名称： PLIN3 Antigen (Recombinant Protein)

别名： perilipin 3; PP17; TIP47; M6PRBP1

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to C terminal 235 amino acids of human PLIN3

### 技术规格

<b>Full name:</b>	perilipin 3
<b>Synonyms:</b>	PP17; TIP47; M6PRBP1
<b>Swissprot:</b>	O60664
<b>Gene Accession:</b>	BC005818
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Mannose 6-phosphate receptors (MPRs) deliver lysosomal hydrolase from the Golgi to endosomes and then return to the Golgi complex. The protein encoded by this gene interacts with the cytoplasmic domains of both cation-independent and cation-dependent MPRs, and is required for endosome-to-Golgi transport.