

兔抗 PAK1(Phospho-Ser199) 多克隆抗体

中文名称：兔抗 PAK1(Phospho-Ser199) 多克隆抗体

英文名称：Anti-PAK1(Phospho-Ser199) rabbit polyclonal antibody

别名：PAKalpha

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：PAK1(Phospho-Ser199)

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The activated kinase acts on a variety of targets. Likely to be the GTPase effector that links the Rho-related GTPases to the JNK MAP kinase pathway. Activated by CDC42 and RAC1. Involved in dissolution of stress fibers and reorganization of focal complexes. Involved in regulation of microtubule biogenesis through phosphorylation of TBCB. Activity is inhibited in cells undergoing apoptosis, potentially due to binding of CDC2L1 and CDC2L2.
Applications:	WB, IHC
Name of antibody:	PAK1(Phospho-Ser199)
Immunogen:	Peptide sequence around phosphorylation site of Serine 19

	9(T-K-S(p)-V-I) derived from Human PAK1.
Full name:	p21 protein (Cdc42/Rac)-activated kinase 1
Synonyms :	PAKalpha
SwissProt:	Q13153
IHC positive control:	Human placenta tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	61 kDa
WB Positive control:	LOVO cells lysates treated with starved
WB Recommended dilution:	500-1000



