

兔抗 PAK1(Phospho-Ser204) 多克隆抗体

中文名称：兔抗 PAK1(Phospho-Ser204) 多克隆抗体

英文名称： Anti-PAK1(Phospho-Ser204) rabbit polyclonal antibody

别名： PAKalpha

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： PAK1(Phospho-Ser204)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a family member of serine/threonine p21-activating kinases, known as PAK proteins. These proteins are critical effectors that link RhoGTPases to cytoskeleton reorganization and nuclear signaling, and they serve as targets for the small GTP binding proteins Cdc42 and Rac. This specific family member regulates cell motility and morphology. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	WB, IHC
Name of antibody:	PAK1(Phospho-Ser204)
Immunogen:	Peptide sequence around phosphorylation site of Serine 20

	4(T-R-S(p)-V-I) derived from Human PAK1.
Full name:	p21 protein (Cdc42/Rac)-activated kinase 1
Synonyms :	PAKalpha
SwissProt:	Q13153
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	61 kDa
WB Positive control:	NIH/3T3 cells lysates treated with UV
WB Recommended dilution:	500-1000



