

兔抗 RUNX1(Phospho-Ser303) 多克隆抗体

中文名称：兔抗 RUNX1(Phospho-Ser303) 多克隆抗体

英文名称： Anti-RUNX1(Phospho-Ser303) rabbit polyclonal antibody

别名： AML1; CBFA2; EVI-1; AMLCR1; PEBP2aB; AML1-EVI-1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： RUNX1(Phospho-Ser303)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

Core binding factor (CBF) is a heterodimeric transcription factor that binds to the core element of many enhancers and promoters. The protein encoded by this gene represents the alpha subunit of CBF and is thought to be involved in the

	<p>development of normal hematopoiesis.</p> <p>Chromosomal translocations involving this gene are well-documented and have been associated with several types of leukemia. Three transcript variants encoding different isoforms have been found for this gene.</p>
Applications:	WB
Name of antibody:	RUNX1(Phospho-Ser303)
Immunogen:	Peptide sequence around phosphorylation site of serine 303 (T-A-F(p)-S-D) derived from Human AML1.
Full name:	runt-related transcription factor 1
Synonyms :	AML1; CBFA2; EVI-1; AMLCR1; PEBP2aB; AML1-EVI-1
SwissProt:	Q01196
WB Predicted band size:	49 kDa
WB Positive control:	Jurkat cells lysates
WB Recommended dilution:	500-1000

