

兔抗 IRF2BP1 多克隆抗体

中文名称：兔抗 IRF2BP1 多克隆抗体

英文名称：Anti-IRF2BP1 rabbit polyclonal antibody

别名：interferon regulatory factor 2 binding protein 1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IRF2BP1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Acts as a transcriptional corepressor in a IRF2-dependent manner; this repression is not mediated by histone deacetylase activities. May act as an E3 ligase towards JDP2, enhancing its polyubiquitination. Represses ATF2-dependent transcriptional activation.
Applications:	ELISA, IHC
Name of antibody:	IRF2BP1
Immunogen:	Synthetic peptide of human IRF2BP1
Full name:	interferon regulatory factor 2 binding protein 1
SwissProt:	Q8IU81
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human esophagus cancer and Human thyroid cancer
IHC Recommend dilution:	30-150



