

兔抗 PAK1 多克隆抗体

中文名称：兔抗 PAK1 多克隆抗体

英文名称：Anti-PAK1 rabbit polyclonal antibody

别名：PAKalpha

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PAK1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a family member of serine/threonine p21-activating kinases, known as PAK proteins. These proteins are critical effectors that link RhoGTPases to cytoskeleton reorganization and nuclear signaling, and they serve as targets for the small GTP binding proteins Cdc42 and Rac. This specific family member regulates cell motility and morphology. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	PAK1

Immunogen:	Synthetic peptide of human PAK1
Full name:	P21 protein (Cdc42/Rac)-activated kinase 1
Synonyms :	PAKalpha
SwissProt:	Q13153
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human liver cancer and Human colon cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	61 kDa
WB Positive control:	293T and K562 cells
WB Recommended dilution:	500-2000



