

## 兔抗 VAV2 多克隆抗体

中文名称：兔抗 VAV2 多克隆抗体

英文名称： Anti-VAV2 rabbit polyclonal antibody

别名： VAV-2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： VAV2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	VAV2 is the second member of the VAV guanine nucleotide exchange factor family of oncogenes. Unlike VAV1, which is expressed exclusively in hematopoietic cells, VAV2 transcripts were found in most tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	VAV2
<b>Immunogen:</b>	Synthetic peptide of humanVAV2
<b>Full name:</b>	vav 2 guanine nucleotide exchange factor

<b>Synonyms :</b>	VAV-2
<b>SwissProt:</b>	P52735
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human breast c ancer
<b>IHC Recommend dilution:</b>	10-50

