

兔抗 IRF2 多克隆抗体

中文名称：兔抗 IRF2 多克隆抗体

英文名称：Anti-IRF2 rabbit polyclonal antibody

别名：interferon regulatory factor 2; IRF-2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：IRF2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	IRF2 encodes interferon regulatory factor 2, a member of the interferon regulatory transcription factor (IRF) family. IRF2 competitively inhibits the IRF1-mediated transcriptional activation of interferons alpha and beta, and presumably other genes that employ IRF1 for transcription activation. However, IRF2 also functions as a transcriptional activator of histone H4.
Applications:	ELISA, WB
Name of antibody:	IRF2
Immunogen:	Synthetic peptide of human IRF2
Full name:	interferon regulatory factor 2
Synonyms:	IRF-2

SwissProt:	P14316
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	39 kDa
WB Positive control:	Hela cell lysate
WB Recommended dilution:	1000-5000

