

## PRELID1 抗原（重组蛋白）

中文名称： PRELID1 抗原（重组蛋白）

英文名称： PRELID1 Antigen (Recombinant Protein)

别名： PRELI domain containing 1; PX19; PRELI; SBBI12; CGI-106

相关类别： 抗原

储存： 冷冻（-20℃）

### 概述

Fusion protein corresponding to a region derived from 36-219 amino acids of human PRELID1

### 技术规格

<b>Full name:</b>	PRELI domain containing 1
<b>Synonyms:</b>	PX19; PRELI; SBBI12; CGI-106
<b>Swissprot:</b>	Q9Y255
<b>Gene Accession:</b>	BC000007
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the late embryogenesis abundant motif-containing protein family. The encoded protein is localized to mitochondria and may function as a cytoprotectant by regulating cell death and differentiation. Alternative splicing results in multiple transcript variants encoding different isoforms. Several related pseudogenes have been identified.

**mlbio** 梅联生物  
Good elisakit producers