

兔抗 RXFP3 多克隆抗体

中文名称：兔抗 RXFP3 多克隆抗体

英文名称：Anti-RXFP3 rabbit polyclonal antibody

别名：relaxin/insulin like family peptide receptor 3; SALPR; RLN3R1; RXFPR3; GPCR135

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：RXFP3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Receptor for RNL3/relaxin-3. Binding of the ligand inhibit cAMP accumulation.Expressed predominantly in brain regions. Highest expression in substantia nigra and pituitary, followed by hippocampus, spinal cord, amygdala, caudate nucleus and corpus callosum, quite low level in cerebellum. In peripheral tissues, relatively high levels in adrenal glands, low levels in pancreas, salivary gland, placenta, mammary gland and testis.
Applications:	ELISA, WB
Name of antibody:	RXFP3

Immunogen:	Synthetic peptide of human RXFP3
Full name:	relaxin/insulin like family peptide receptor 3
Synonyms:	SALPR; RLN3R1; RXFPR3; GPCR135
SwissProt:	Q9NSD7
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	51 kDa
WB Positive control:	K562, 293T and NIH/3T3 cell lysates
WB Recommended dilution:	200-1000

