

兔抗 RXFP4 多克隆抗体

中文名称: 兔抗 RXFP4 多克隆抗体

英文名称: Anti-RXFP4 rabbit polyclonal antibody

别名: relaxin family peptide/INSL5 receptor 4; GPR100; RLN3R2; RXFP4; GPCR142

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: RXFP4

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	High affinity receptor for INSL5. Also acts as receptor for relaxin or RLN3/relaxin-3, as well as bradykinin and kallidin. Binding of the ligand inhibits cAMP accumulation.
Applications:	ELISA, IHC
Name of antibody:	RXFP4
Immunogen:	Synthetic peptide of human RXFP4
Full name:	relaxin family peptide/INSL5 receptor 4
Synonyms:	GPR100; RLN3R2; RXFP4; GPCR142
SwissProt:	Q8TDU9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer

IHC Recommend dilution:	50-100
-------------------------	--------

