

兔抗 PAK2(Phospho-Ser197) 多克隆抗体

中文名称：兔抗 PAK2(Phospho-Ser197) 多克隆抗体

英文名称： Anti-PAK2(Phospho-Ser197) rabbit polyclonal antibody

别名： PAK65; PAKgamma

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： PAK2(Phospho-Ser197)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The p21 activated kinases (PAK) are critical effectors that link Rho GTPases to cytoskeleton reorganization and nuclear signaling. The PAK proteins are a family of serine/threonine kinases that serve as targets for the small GTP binding proteins, CDC42 and RAC1, and have been implicated in a wide range of biological activities. The protein encoded by this gene is activated by proteolytic cleavage during caspase-mediated apoptosis, and may play a role in regulating the apoptotic events in the dying cell.
Applications:	IHC
Name of antibody:	PAK2(Phospho-Ser197)

Immunogen:	Peptide sequence around phosphorylation site of Serine 197(T-R-S(p)-V-I) derived from Human PAK2.
Full name:	p21 protein (Cdc42/Rac)-activated kinase 2
Synonyms :	PAK65; PAKgamma
SwissProt:	Q13177
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100

