

## 兔抗 VCAM1 多克隆抗体

中文名称: 兔抗 VCAM1 多克隆抗体

英文名称: Anti-VCAM1 rabbit polyclonal antibody

别名: CD106; INCAM-100

抗原: VCAM1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the Ig superfamily and encodes a cell surface sialoglycoprotein expressed by cytokine-activated endothelium. This type I membrane protein mediates leukocyte-endothelial cell adhesion and signal transduction, and may play a role in the development of atherosclerosis and rheumatoid arthritis. Three alternatively spliced transcripts encoding different isoforms have been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	VCAM1
<b>Immunogen:</b>	Fusion protein of human VCAM1
<b>Full name:</b>	vascular cell adhesion molecule 1

<b>Synonyms :</b>	CD106; INCAM-100
<b>SwissProt:</b>	P19320
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human prostate cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

