

## 兔抗 IRF6 多克隆抗体

中文名称：兔抗 IRF6 多克隆抗体

英文名称：Anti-IRF6 rabbit polyclonal antibody

别名：LPS; PIT; PPS; VWS; OFC6; PPS1; VWS1

抗原：IRF6

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the interferon regulatory transcription factor (IRF) family. Family members share a highly-conserved N-terminal helix-turn-helix DNA-binding domain and a less conserved C-terminal protein-binding domain. The encoded protein may be a transcriptional activator. Mutations in this gene can cause van der Woude syndrome and popliteal pterygium syndrome. Mutations in this gene are also associated with non-syndromic orofacial cleft type 6. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	IRF6

<b>Immunogen:</b>	Fusion protein of human IRF6
<b>Full name:</b>	interferon regulatory factor 6
<b>Synonyms :</b>	LPS; PIT; PPS; VWS; OFC6; PPS1; VWS1
<b>SwissProt:</b>	O14896
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human cervical cancer and human lung cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	53 kDa
<b>WB Positive control:</b>	Mouse kidney and lung tissue
<b>WB Recommended dilution:</b>	200-1000



