

## 兔抗 VCAN 多克隆抗体

中文名称: 兔抗 VCAN 多克隆抗体

英文名称: Anti-VCAN rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

抗 原: VCAN

宿 主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the aggrecan/versican proteoglycan family. The protein encoded is a large chondroitin sulfate proteoglycan and is a major component of the extracellular matrix. This protein is involved in cell adhesion, proliferation, proliferation, migration and angiogenesis and plays a central role in tissue morphogenesis and maintenance. Mutations in this gene are the cause of Wagner syndrome type 1. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	VCAN
<b>Immunogen:</b>	Fusion protein of human VCAN

<b>Full name:</b>	versican
<b>Synonyms :</b>	WGN; ERVR; GHAP; PG-M; WGN1; CSPG2
<b>SwissProt:</b>	P13611
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	10-50

