

兔抗 S1PR3 多克隆抗体

中文名称：兔抗 S1PR3 多克隆抗体

英文名称：Anti-S1PR3 rabbit polyclonal antibody

别名：EDG3; LPB3; S1P3; EDG-3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：S1PR3

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the EDG family of receptors, which are G protein-coupled receptors. This protein has been identified as a functional receptor for sphingosine 1-phosphate and likely contributes to the regulation of angiogenesis and vascular endothelial cell function. Receptor for the lysosphingolipid sphingosine 1-phosphate (S1P). S1P is a bioactive lysophospholipid that elicits diverse physiological effects on most types of cells and tissues. When expressed in rat HTC4 hepatoma cells, is capable of mediating S1P-induced cell proliferation and suppression of

	f apoptosis.
Applications:	ELISA, WB, IHC
Name of antibody:	S1PR3
Immunogen:	Synthetic peptide of human S1PR3
Full name:	Sphingosine-1-phosphate receptor 3
Synonyms :	EDG3; LPB3; S1P3; EDG-3
SwissProt:	Q99500
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human cervical cancer and Human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	42 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	1000-5000



