

兔抗 S100A6 多克隆抗体

中文名称：兔抗 S100A6 多克隆抗体

英文名称： Anti-S100A6 rabbit polyclonal antibody

别名： 2A9; PRA; 5B10; CABP; CACY

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： S100A6

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the S100 family of proteins containing 2 EF-hand calcium-binding motifs. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells, and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100 genes include at least 13 members which are located as a cluster on chromosome 1q21.
Applications:	ELISA, WB
Name of antibody:	S100A6

Immunogen:	Fusion protein of human S100A6
Full name:	S100 calcium binding protein A6
Synonyms :	2A9; PRA; 5B10; CABP; CACY
SwissProt:	P06703
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	10 kDa
WB Positive control:	SKOV3 and A431 cells
WB Recommended dilution:	500-2000

