

## 兔抗 pan CDH 多克隆抗体

- 中文名称：兔抗 pan CDH 多克隆抗体
- 英文名称：Anti-pan CDH rabbit polyclonal antibody
- 相关类别：一抗
- 储 存：冷冻（-20℃）
- 宿 主：Rabbit
- 抗 原：pan CDH
- 反应种属：Human, Mouse, Rat
- 标 记 物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a classical cadherin from the cadherin superfamily. The encoded protein is a calcium dependent cell-cell adhesion glycoprotein comprised of five extracellular cadherin repeats, a transmembrane region and a highly conserved cytoplasmic tail. Mutations in this gene are correlated with gastric, breast, colorectal, thyroid and ovarian cancer. Loss of function is thought to contribute to progression in cancer by increasing proliferation, invasion, and/or metastasis. The ectodomain of this protein mediates bacterial adhesion to mammalian cells and the cytoplasmic domain is required for internalization. Identified transcript variants arise from mutation at consensus splice sites.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	pan CDH
<b>Immunogen:</b>	Synthetic peptide of human pan CDH
<b>Full name:</b>	pan-cadherin
<b>SwissProt:</b>	P12830
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human cervical cancer and Human lung cancer
<b>IHC Recommend dilution:</b>	25-100

