

## 兔抗 S100A14 多克隆抗体

中文名称：兔抗 S100A14 多克隆抗体

英文名称： Anti-S100A14 rabbit ployclonal antibody

别名： BCMP84; S100A15

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： S100A14

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the S100 protein family which contains an EF-hand motif and binds calcium. The gene is located in a cluster of S100 genes on chromosome 1. Levels of the encoded protein have been found to be lower in cancerous tissue and associated with metastasis suggesting a tumor suppressor function
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	S100A14
<b>Immunogen:</b>	Fusion protein of human S100A14
<b>Full name:</b>	S100 calcium binding protein A14

<b>Synonyms :</b>	BCMP84; S100A15
<b>SwissProt:</b>	Q9HCY8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human liver cancer and human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100

