

兔抗 IRS1 (Phospho-Ser312)多克隆抗体

中文名称： 兔抗 IRS1 (Phospho-Ser312)多克隆抗体

英文名称： Anti-IRS1 (Phospho-Ser312) rabbit polyclonal antibody

别名： HIRS-1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： IRS1 (Phospho-Ser312)

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	May mediate the control of various cellular processes by insulin. When phosphorylated by the insulin receptor binds specifically to various cellular proteins containing SH2 domains such as phosphatidylinositol 3-kinase p85 subunit or GRB2. Activates phosphatidylinositol 3-kinase when bound to the regulatory p85 subunit
Applications:	WB
Name of antibody:	IRS1 (Phospho-Ser312)
Immunogen:	Synthetic peptide of human IRS1 (Phospho-Ser312)
Full name:	insulin receptor substrate 1 (Phospho-Ser312)
Synonyms :	HIRS-1

SwissProt:	P35568
WB Predicted band size:	180 kDa
WB Positive control:	293 cells treated with insulin
WB Recommended dilution:	500-1000

