

PTPN3 抗原(重组蛋白)

中文名称: PTPN3 抗原(重组蛋白)

英文名称: PTPN3 Antigen (Recombinant Protein)

储 存: 冷冻(-20℃)

别 名: protein tyrosine phosphatase, non-receptor type 3; PTPH1; PTP-H1

相关类别: 抗原

概述

Fusion protein corresponding to a region derived from 29-312 amino acids of human PTPN3

技术规格

Full name:	protein tyrosine phosphatase, non-receptor type 3
Synonyms:	PTPH1; PTP-H1
Swissprot:	P26045
Gene Accession:	BC126117
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a member of the protein tyr osine phosphatase (PTP) family. PTPs are known to be signaling molecules that regulate a variety of cellular processes including ce II growth, differentiation, mitotic cycle, and oncogenic transformati on. This protein contains a C-terminal PTP domain and an N-term inal domain homologous to the band 4.1 superfamily of cytoskele tal-associated proteins.