

## 兔抗 PARD6A 多克隆抗体

中文名称：兔抗 PARD6A 多克隆抗体

英文名称： Anti-PARD6A rabbit polyclonal antibody

别名： PAR6; PAR6C; TAX40; PAR-6A; TIP-40; PAR6alpha

抗原： PARD6A

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the PAR6 family and encodes a protein with a PSD95/Discs-large/ZO1 (PDZ) domain and a semi-Cdc42/Rac interactive binding (CRIB) domain. This cell membrane protein is involved in asymmetrical cell division and cell polarization processes as a member of a multi-protein complex. The protein also has a role in the epithelial-to-mesenchymal transition (EMT) that characterizes the invasive phenotype associated with metastatic carcinomas. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	PARD6A
<b>Immunogen:</b>	Fusion protein of human PARD6A
<b>Full name:</b>	par-6 family cell polarity regulator alpha
<b>Synonyms :</b>	PAR6; PAR6C; TAX40; PAR-6A; TIP-40; PAR6alpha
<b>SwissProt:</b>	Q9NPB6
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human colon cancer and human breast cancer
<b>IHC Recommend dilution:</b>	50-200

