

兔抗 VEGFB 多克隆抗体

中文名称：兔抗 VEGFB 多克隆抗体

英文名称：Anti-VEGFB rabbit polyclonal antibody

别名：VRF, VEGFL

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：VEGFB

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the PDGF (platelet-derived growth factor)/VEGF (vascular endothelial growth factor) family. The VEGF family members regulate the formation of blood vessels and are involved in endothelial cell physiology. This member is a ligand for VEGFR-1 (vascular endothelial growth factor receptor 1) and NRP-1 (neuropilin-1). Studies in mice showed that this gene was co-expressed with nuclear-encoded mitochondrial genes and the encoded protein specifically controlled endothelial uptake of fatty acids. Alternatively spliced transcript variants encoding distinct isoforms have been identified.

Applications:	ELISA, WB, IHC
Name of antibody:	VEGFB
Immunogen:	Synthetic peptide of human VEGFB
Full name:	vascular endothelial growth factor B
Synonyms :	VRF, VEGFL
SwissProt:	P49765
ELISA Recommended dilution:	5000-240000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	21 kDa
WB Positive control:	Mouse brain tissue
WB Recommended dilution:	500-2000



