

兔抗 SAE1 多克隆抗体

中文名称：兔抗 SAE1 多克隆抗体

英文名称： Anti-SAE1 rabbit polyclonal antibody

别名： AOS1; SUA1; UBLE1A; HSPC140

相关类别： 一抗

储存： 冷冻 (-20℃)

抗原： SAE1

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Posttranslational modification of proteins by the addition of the small protein SUMO, or sumoylation, regulates protein structure and intracellular localization. SAE1 and UBA2 form a heterodimer that functions as a SUMO-activating enzyme for the sumoylation of proteins.
Applications:	ELISA, WB, IHC
Name of antibody:	SAE1
Immunogen:	Fusion protein of human SAE1
Full name:	SUMO1 activating enzyme subunit 1
Synonyms :	AOS1; SUA1; UBLE1A; HSPC140

SwissProt:	Q9UBE0
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human cervical cancer
IHC Recommend dilution:	15-50
WB Predicted band size:	38 kDa
WB Positive control:	293T and A431 cells
WB Recommended dilution:	500-2000



