

## 兔抗 IRX1 多克隆抗体

中文名称: 兔抗 IRX1 多克隆抗体

英文名称: Anti-IRX1 rabbit polyclonal antibody

别名: iroquois homeobox 1; IRX-5; IRXA1

抗原: IRX1

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗


标记物: Unconjugate


克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the Iroquois homeobox protein family. Homeobox genes in this family are involved in pattern formation in the embryo. The gene product has been identified as a tumor suppressor in gastric and head and neck cancers. A pseudogene of this gene is located on chromosome 13.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	IRX1
<b>Immunogen:</b>	Synthetic peptide of human IRX1
<b>Full name:</b>	iroquois homeobox 1
<b>Synonyms:</b>	IRX-5; IRXA1

<b>SwissProt:</b>	P78414
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	40-200
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	Human heart tissue lysate
<b>WB Recommended dilution:</b>	500-2000





**mlbio** 梅联生物  
Good elisakit producers

