

兔抗 SAMD7 多克隆抗体

中文名称：兔抗 SAMD7 多克隆抗体

英文名称： Anti-SAMD7 rabbit ployclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SAMD7

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	SAMD7 (Sterile alpha motif domain-containing protein 7), is a 446 amino acid protein encoded by the SAMD7 gene which maps to human chromosome 3. Chromosome 3 is made up of about 214 million bases encoding over 1,100 genes, including a chemokine receptor (CKR) gene cluster and a variety of human cancer-related gene loci. Key tumor suppressing genes on chromosome 3 include those that encode the apoptosis mediator RASSF1, the cell migration regulator HYAL1 and the angiogenesis suppressor SEMA3B.
Applications:	ELISA, IHC
Name of antibody:	SAMD7

Immunogen:	Fusion protein of human SAMD7
Full name:	sterile alpha motif domain containing 7
SwissProt:	Q7Z3H4
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human breast cancer
IHC Recommend dilution:	25-100

