

RAB20 抗原（重组蛋白）

中文名称： RAB20 抗原（重组蛋白）

英文名称： RAB20 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

别 名： RAB20, member RAS oncogene family

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-234 amino acids of human RAB20

技术规格

Full name:	RAB20, member RAS oncogene family
Swissprot:	Q9NX57
Gene Accession:	BC026025
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Plays a role in apical endocytosis/recycling. Plays a role in the maturation and acidification of phagosomes that engulf pathogens, such as S.aureus and M.tuberculosis. Plays a role in the fusion of phagosomes with lysosomes.