

兔抗 SAP30BP 多克隆抗体

中文名称：兔抗 SAP30BP 多克隆抗体

英文名称： Anti-SAP30BP rabbit polyclonal antibody

别名： SAP30 binding protein; HTRG; HTRP; HCNGP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SAP30BP

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Induces cell death. May act as a transcriptional corepressor or of a gene related to cell survival. May be involved in the regulation of beta-2-microglobulin genes.
Applications:	ELISA, WB, IHC
Name of antibody:	SAP30BP
Immunogen:	Fusion protein of human SAP30BP
Full name:	SAP30 binding protein
Synonyms:	HTRG; HTRP; HCNGP
SwissProt:	Q9UHR5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer

IHC Recommend dilution:	50-200
WB Predicted band size:	34 kDa
WB Positive control:	HUVEC cell lysate
WB Recommended dilution:	500-2000



